miccra

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Labortechnik | Laboratory technology | 实验室仪器

Technikum | Pilot Plants | 中试产品

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The „ART“ of homogenizing

As a specialist in the homogenizing technology, ART Prozess- & Labortechnik GmbH & Co. KG is developing and manufacturing powerful, high quality and high speed dispersing- and homogenizing instruments since last Century.

Our ambition has always been to offer easy to handle instruments for achieving reliable, reproducible sample results coupled with an unmatched performance to price ratio.

The MICCRA-Series of ART disperses, homogenizes, emulsifies, suspends, cuts up and disrupts sample materials of nearly all kinds, based on the rotor-stator principle, but with higher performance.

The modular construction of the MICCRA-Series allows flexible and multiple combinations of different powerful drives, shafts in several lengths, and all the dispersing accessories, tailored to the individual requirements.

The systems can either be operated as a hand-held, post or wall mounted or a tripod support device. With our goal of ensuring customer satisfaction, we emphasize our own R & D and the cooperation with leading universities and research institutes.

Continually in touch with our customers, we develop and put into practice solutions matching the specific requirements of users.

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By a light and compact design, all lab-drives of the MICCRA-Series can be used as hand-held units.

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High-Speed-drives with continuously adjustable rpm-rates up to 39,000 rpm, reduce processing times significantly without thermal stressing of the sample material.

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High power for handling large volumes and high viscosities. Unique application ranges from 0.1 ml up to 100l.

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Comfortable and safe working with smooth start of the drive.

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
Uniform shaft size:
Every drive fits to every tool combination of ART.

主机
The drive

刀头
The tool

套管 S 20
Shaft tube S 20

轴 W 20 含 PTFE
轴承 P20
Axle W 20 with
PTFE-bearing
P 20

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High precision tools for highest reliable results at highest peripheral speeds.

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Due to a uniform shaft size, every tool combination fits to every drive within the MICCRA-series.

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Flexible modular system for high economy and small subsequent costs. By combination of different rotors and stators the tool can be adapted to nearly each range of application w.o. large expenditure.

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Quick-Change-System: All elements, rotors and stators can be exchanged with a flick of the wrist in no time at all.

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All tools of the MICCRA-Series are made of high-quality AISI 316L stainless steel. Easy to clean and long-lasting.

Upon request our tools are also available in titanium, Hastelloy, Inconel as well as other materials.

MICCRA

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Dispersing and homogenizing drives

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The economy drive MICCRA D-1 is a small hand-held dispersing and homogenizing device for volumes from 0.1 up to 500 ml. It was especially developed for volumes treated in the bio- and gene technology.

It's an economic drive for a small price with a „soft-touch“ housing. More powerful and high-speed than comparable systems on the market.

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Light, compact, strong and extremely high speed.

These are the topics of this extraordinary drive which fits to all common laboratory applications from the agrarian technology to cell biology.

A unique volume range from 0.1 up to 20.000 ml can be handled with this unit and its unequalled results.

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Despite its small size, this drive is a real high power unit.

The successor of the MICCRA D-13, with its integrated constant electronics and its more powerful 1520 Watt motor works out large volumes, keeping a constant speed even with increasing viscosities.

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Volume range (depending on vessel to be used)

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Our „PICO“ biotech tool range applies particularly for smaller sample quantities in the bio and gene-technological sample preparation, especially in PCR diagnostics. Whether you would like to extract RNA or cut up fabric for other analysis procedures.

The tools of the PICO biotech- series will accelerate and simplify your work without large sample heating up.

They are also reliable aids for cell explanations and enzyme studies. Combined with the high numbers of revolutions of the drives MICCRA D-1 and D-9 you will receive excellent results within seconds, which have been confirmed to us by numerous laboratories.

刀头类型 Tool type

定子类型 Stator type

定子外径 Outer diameter of stator

转子类型 Rotor type

转子外径 Outer diameter of rotor

最大剪切速度 Max. peripheral speed

刀头长度 Tool length

最小浸入深度 Max./min. insertion depth

处理量 (水) Treatable Volume (acc. to H2O)

DS-5/K 1

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Dispersing tool for small volumes

刺钉刀头 spike head

5 mm

剪切转子 Cutting rotor

3,6 mm

7,3 m/s

115 mm

30 mm / 3 mm

0,1 - 20 ml

DS-8/P

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Dispersing tool for small volumes

刺钉刀头 spike head

8 mm

剪切转子 Cutting rotor

5,6 mm

11,4 m/s

155 mm

80 mm / 20 mm

1 - 100 ml

DS-14/P

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Dispersing tool for small volumes

刺钉刀头 spike head

14 mm

剪切转子 Cutting rotor

10,6 mm

21,6 m/s

200 mm

125 mm / 30 mm

10 - 500 ml

Vario

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Dispersing and homogenizing tools - BATCH mode

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The application range in science, research and production is various. Due to our „quick-change-system“, the disassembly and change are very easy.

The modular system results in cost savings and frees up working space. Due to newly developed rotors even better results are achieved, even by a cutback of our tool diversity.

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
2 in 1 tool with cutter for fibrous materials

定子类型 Stator type

粗的 coarse

精细 fine

精细 fine

刺钉刀头 spike head

定子外径 Outer diameter of stator

20 mm

20 mm

20 mm

20 mm

转子类型 Rotor type

剪切-混合-转子 Cutting-Mixing-rotor

剪切-混合-转子 Cutting-Mixing-rotor

乳化转子 Emulsifying rotor

剪切-混合-转子 Cutting-Mixing-rotor

转子外径 Outer diameter of rotor

15 mm

15 mm

15 mm

15 mm

最大剪切速度 Max. peripheral speed

30,6 m/s

30,6 m/s

30,6 m/s

30,6 m/s

刀头长度 Tool length

250 mm

250mm

250mm

250 mm

最小浸入深度 Max./min. insertion depth

200 mm / 40 mm

200mm / 40mm

200mm / 40mm

200 mm / 40 mm

处理量（水） Treatable Volume (acc. to H2O)

20 - 5.000 ml

20 - 3.000ml

20 - 3.000ml

20 - 3.000 ml

Vario

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Dispersing and homogenizing tools - BATCH mode

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Highly polished standard AISI 316L SS dispersing tools ensure a long life . All elements can be assembled and disassembled quickly, and the excellent surface quality guarantees you a simple cleaning of the components.

All tools are also available in XL (tool length 327 mms) and XXL version (420 mms).

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Our economic „2/P-bundles“ consist of a coarse and a fine stator each, as well as a rotor.

The desired end fineness is achieved by changing the stator from coarse to fine without large expenditure by the „quick-change-system“.

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All-in-one tool for pasty, fibrous and liquid components

All-in-one tool for pasty, fibrous and liquid components

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
Get the latest product info
directly on your smartphone!

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

.....
With our VARIOMax - line for large volumes, you are able to mix volumes up to 100 l easily.

Our standard lab-drive MICCRA D-15 provides, on this occasion, constant speeds even with increasing viscosity.

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
Dispersing and homogenizing system - InLine mode

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The flexible solution for short-term operation. Our standard laboratory homogenizers MICCRA D-15 or MICCRA D-9 provide the needed drive for your work.

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The flow chamber system made for continuous operation in lab and pilot plant scale. Firmly connected with our MICCRA D-27 high-end drive, you will achieve unique product results up to nanometer range.

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Safe, flexible, efficient.

The MICCRA DFK Lab is an easy and quick to clean, efficient flow chamber system for your lab and small scale production purposes.

As a result of the modular construction, all known rotor/stator-components from the batch-series are uncompromisingly available for this chamber also.

Due to the significant higher circumferential speeds, end-finenesses up to 40 nm can be achieved according to product features and application.

As yet, this was not possible with such small systems.

推荐主机 | Recommended drives

Miccra D-9 短时间工作 | (short-time operation)

Miccra D-15 短时间工作 | (short-time operation)

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
The DFK system allows you the processing of larger volumes in InLine-/
flowthrough-procedure in lab and pilot plant scale.

The solid construction method made of high-quality stainless steel and the highly
resistant sealing materials guarantee a versatile application in all industrial
branches to:

- chemically and/or physically change of media capable of flow
- carry out smaller and therefore cost-saving attempts in the
field of research
- create e.g. new recipes,
- transfer the gained insights in the lab scale to your future
production sizes (scale up)

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2 sealing types

You have the choice between a single
mechanical seal and a lip seal with FDA
approval.

The patent pending cartridge construction
provides easy cleaning of the lip seal and
may even be replaced by the user itself, if
necessary. Therefore, you save costs and
time for repair work.

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Flexible use of different rotor / stator
configurations and sizes

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Revolutionary sealing system in cartridge
construction allows quick and easy repla-
cement of the sealing lips

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Standard high-polished tool components

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Accessories such as stainless steel containers with different volumes
(1, 3 or 5 liters) in single or double-walled construction, hose nipples,
fittings and three way valves etc. are standardly available. Further
volumes/sizes available upon request.

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The ideal solution for short-time operation

The DFK 1.0 flow chamber system is the entry to a fast and cost-effective product development. Again, you have the possibility of being able to use the batch series of our known modular rotor / stator components without compromise

Its usefulness as e.g. a product developer lies in the following facts..

- ... That for the first time, you can fall back on a system which is absolutely the same structure as it is in the batch operation of the MICCRA laboratory series.
- ... That you can flexibly attach pipes, fittings, tanks (whether single or double-walled, closed, open), valves of any kind via Tri-Camp connections to the system and simulate your production more or less on a small scale.
- ... That you even can achieve product endfinenesses with a mechanical system, which lie in the nanometer range.

And at a cost that is unmatched. The gained knowledge can then be adopted into the next larger scale in pilot plant or production.

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Product optimization

The inline flow chamber system DFK 1.4 is the high-end solution for continuous operation in research, development and small production scale. Based on the MICCRA D-27, the powerful 2.7 kW induction motor with speeds up to 36,000 1/min, product results with a much higher quality and efficiency are possible.

The results are reproducible at any time and a scale-up to the next higher volume areas is possible.

To mixing and size reduction effects is to optimize a wide variety of rotor / stator combinations.

Thanks to a bearingless drive shaft between the flow chamber and motor drive the system is characterized by its remarkable smoothness, and is particularly easy to clean and maintain.

Thus, the ideal combination between MICCRA D-27 and the DFK 1.4 By flow chamber system as an ideal equipment for their crossover products.

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The new high-end homogenizer, built for the highest requirements in research and development.

Its powerful 2.7 kW stainless steel three-phase asynchronous motor with up to 36,000 rpm delivers product results up to nanoscale, at a surprisingly low noise rate (63 dB(A)).

The patent-pending, unique cartridge sealing system (in connection with the DFK-flow chamber) withstands also pressure and vacuum applications easily and is particularly hassle-free to maintain.*

The perfect solution for the continuous operation in laboratory and pilot plant, hands down.

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Motor type:
Voltage:
Electricity:
Torque:
Power:
Frequency:
Speed regulation:
Speed range:
Housing type:
Drive type:
Safety type:
Cooling type:
Sealing air:
Sealing type:
Dimensions:
Weight:

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three phase asynchronous
185V 3~
max. 10,6 A
max. 200 Ncm
max. 2700 W
50 - 400 Hz (600 Hz)
via frequency converter
3.000 - 24.000 rpm. (optional up to 36.000 rpm)
Stainless steel
S1 (continuous operation)
IP 54
integrated water cooling circuit
integrated compressed air supply
lubricated-for-life spindle ball bearing
112 x 112 x 442mm (WxHxD)
9.5kg

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HFU-27 high frequency converter

High frequency converter for MICCRA D-27 asynchronous drive

Compact and resistant table and / or mounted device with integrated digital display and short circuit monitoring. Use in the public mainsupply by integrated PFC (Power Factor Correction).

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Input (AC):
Output (AC):
Power output:
Modulation:
Nominal current:
Peak current:
Safety type:
Ambient temp.:
Equipment:
Operation type
Dimensions:
Weight:

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1 x 115 to 1x230V 50/60Hz
3x 0 to 220V
3,8kVA (3,8kW)
Pulse width modulation
10A
28A
IP 20
5 to 40°C
delivery with connector set
S1 (continuous operation)
120x275x220mm (HxWxD)
4.9 kgs

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Motor connection:
USB connection:
Power supply:

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6-pin lead-out terminal
yes
On-/off switch on front
LCD status display for
displaying operational
data such as:
parameters, speed
setpoint / actual, load
indicator, active drive,
error messages

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
LED status indicator:
Ready, active load detector,
keypad for operation,
drive selection and
parameterization

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for fixing of vessels up to diameter 80 mm,
incl. boss head KM20.

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for fixing of cylindrical vessels
incl. boss head KM20.

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For stationary work
stainless-steel,
surface dimensions: 300x210mm
incl. 2 rods

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for stationary work
stainless steel- construction

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Tool fitting, resp. holding device,
PTFE for dispersing with diameter 20 mm,
with adapter for ground taper socket NS 29

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for adapting of tool fitting WA20
on reactor lids such as BÜCHI-
reactors

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for fixing of accessories and units on
plate stand, heavy duty version

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for dispersing tools
OW1 for EMR20, 30, 40
OW2 for SMIR 20
OW3 for SMIR 30, 40

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for a precise adjustment of process
time down to 1 sec.
Up to 4 devices connectable

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Special heavy duty glassware for a
quick, intense and safe sample pre-
paration. Different sizes available

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Please also ask for our separate brochures for the following ranges:

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How special your application may be, we might have the right solution for you.

Just ask our skilled professionals for custom solutions made by ART.

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Handle volumes of 50-500 l easily in batch, and flow rates up to 10t / h inline-operation with our „Micra Pi“ pilot plant series.


The seamless transition from lab to production scale.

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Our „Micra Pro“ series offers professional and custom-made solutions in terms of production scale mixing and dispersing.

安全箱
Safety cabin


批次式中试设备
Batch operated pilot plant

在线式中试设备
Inline operated pilot plant

批次式/在线式工业设备
Inline/batch operated production plant

miccra